

Overview

The assessment of a patient`s comprehensive lipoprotein profile gives insights in cardiovascular health and risk for e.g. diabetes, metabolic syndrom or stroke.

Regulatory Status

RESEARCH USE ONLY (RUO): The RUO-Test contains no medical diagnosis and does not replace the medical diagnosis of a physician.

Method

lifespın LipoPRO V 2.0.0 based on nuclear magnetic resonance (NMR) data.

Test Parameter

Main Lipid Parameters											
• Total Cholesterol • Total Triglycerides		• LDL-Cholesterol • HDL-Cholesterol		• Apo-A1 / A2 • Apo-B100							
Particle Concentrations		Particle Sizes		Parameter Ratios							
• Lipoprotein Main <u>Fraction</u> • Non-HDL <u>particle number</u> • LDL-p • IDL-p • VLDL-p • HDL-p • Chylomicrons-p • HDL-Subfraction • HDL-1-p to HDL-4-p • LDL-Subfraction • LDL-1-p to LDL-6-p		• Mean <u>diameter of Chylomicrons</u> • Mean diameter of VLDL particles • Mean <u>diameter of VLDL-1 to VLDL-5 particles</u> • Mean <u>diameter of LDL particles</u> • Mean <u>diameter of LDL-1 to LDL-6 particles</u> • Mean <u>diameter of HDL particles</u> • Mean <u>diameter of HDL-1 to HDL-4 particles</u>		• LDL-Cholesterol/ HDL-Cholesterol • Apo-B100 / Apo-A • Free LDL-Cholesterol / LDL-Cholesterol • Free HDL-Cholesterol / HDL-Cholesterol							
Concentrations in Main & Subfractions											
• Cholesterol in • HDL • LDL • IDL • VLDL • Chylomicrons • LDL-1 to LDL-6 • HDL-1 to HDL-4 • VLDL-1 to VLDL-5		• Free Cholesterol in • HDL • LDL • IDL • VLDL • Chylomicrons • LDL-1 to LDL-6 • HDL-1 to HDL-4 • VLDL-1 to VLDL-5		• Esterified Cholesterol in • HDL • LDL • IDL • VLDL • Chylomicrons • LDL-1 to LDL-6 • HDL-1 to HDL-4 • VLDL-1 to VLDL-5		• Triglycerides in • HDL • LDL • IDL • VLDL • Chylomicrons • LDL-1 to LDL-6 • HDL-1 to HDL-4 • VLDL-1 to VLDL-5		• Phospholipids in • HDL • LDL • IDL • VLDL • Chylomicrons • LDL-1 to LDL-6 • HDL-1 to HDL-4 • VLDL-1 to VLDL-5		• Apo-A1/A2 in • HDL • Apo-B100 in • LDL • IDL • VLDL • Chylomicrons	

Specimen Type

Matrix:

Human blood serum.

Specimen Requirements

Patient Preparation:

Before blood sampling patients should be fasting for 8 hours.

Collection Container / Tube:

Sarstedt S-Monovette (red top). Serum gel barrier tubes are not acceptable.

Specimen Volume: 1 mL.

Specimen Minimum Volume: 0.5 mL.

Collection instruction:

- Allow isopropyl alcohol (from phlebotomy site prep) to dry thoroughly before venipuncture.
- Centrifuge and aliquot serum into an appropriate plastic vial.

Storage Instruction:

Send-in-Services:

Keep frozen and temperature-guided shipment on dry ice.

Reject Due To

- Specimen other than serum
- Improper labeling
- Serum collected using gel barrier tube
- Specimen received in inappropriate container
- Gross hemolysis
- Gross lipemia
- Gross icterus
- Precipitation

Method Description

The analysis software “lifespın LipoPRO V1.0.0.cloud” is used for the quantification of lipoprotein particles in serum based on NMR data.

Limitations

In case of pathological Triglyceride values, the precision of LDL-Cholesterol-value might be effected.

The parameter values can also be impaired if the blood is heavily contaminated with non- or-naturally occurring substances.

Expected Turnaround Time

Send-in-Services

The processing of the order begins at the earliest with the arrival of the sample material relevant for the analyses and the associated information. Time of analyses, capacities and delivery times by arrangement (in General: up to 200 Samples in 3 working days).

Performing Laboratory Location
Regensburg, Germany

Software as a Service (SaaS)

Human samples measured on an inhouse NMR spectrometer qualified by lifespın can be processed via cloud solution. (turnaround time is 500 samples in less than 25 minutes).

Performing Laboratory Location
Customer Site

Useful for

Overall assessment of the current status of cardiovascular health

According to the World Health Organization, cardiovascular disease is one of the most common causes of death worldwide. Approximately one third of all deaths are caused by cardiovascular disease.

Evaluation of lipoprotein particles has been used to support management of cardiovascular disease (CVD) risk for over 15 years, and lipoprotein subclass analysis has become a valuable tool to help clinicians better stratify patients at risk.

The lifespın LipoPRO-test is used to help assess the risk for CVD in patients based on traditional or emerging risk factors such as quantification of the number and size of lipoprotein particles.

If both standard and emerging lipid parameters such as number and size rated green, no hint of coronary artery disease is indicated.

If the parameters are rated orange or red, closer or further considerations are recommended.

According to peer reviewed literature lipoprotein subclasses hint to the following health conditions:

- When LDL-Cholesterol is normal and LDL-P is elevated (lipid discordance) CVD risk tracks with levels of LDL-P. LDL-P is more strongly associated with risk for CVD events and atherosclerosis than LDL-C.
- LDL-P can identify patients with intermediate risk for CVD events. For intermediate-risk patients, LDL-P provides additional insight to CVD risk over above standard risk factors.
- HDL-P is inversely associated with CVD risk. HDL-P is associated with reduced CVD risk, with or without statin therapy.
- Lipid fraction, including small LDL-P, LDL Size, Large HDL-P, HDL Size, Large VLDL-P, and VLDL Size, are correlated with risks of metabolomic syndrome, diabetes, coronary artery disease and stroke.

Relevant Ranges

Parameter	Reference Ranges [mg/dL]	Result
Cholesterol	<200	Desirable
	200-239	Borderline High
	>239	High
LDL-Cholesterol	<70	Best for people who have coronary artery disease.
	<100	Optimal for people at risk of coronary artery disease or who have diabetes.
	100-129	Near Optimal if there is no coronary artery disease. High if there is coronary disease.
	129-159	Borderline High if there is no coronary artery disease. High if there is coronary disease.
	159-189	High if there is no coronary artery disease. Very high if there is coronary disease.
	>189	Very High
HDL-Cholesterol	<40	Poor
	40-59	Better
	>59	Best
Triglyceride	<150	Desirable
	150-199	Borderline High
	199-499	High
	>499	Very High

The used relevant ranges are sexspecific clinical references from the guideline of Mayo Clinic, USA.

Reference Ranges*

Parameter	Reference Ranges (f) mg/dL	Reference Ranges (m) mg/dL	Reference Ranges (a) mg/dL
Apo-A1	147.2-269.3	125.7-240.8	130.4-262.1
Apo-A2	31.4-54.5	27.1-49.7	29.1-53.1
Apo-B100	48.7-162.2	47.7-179.1	48.1-170.6

Ratios

Parameter	Reference Ranges (f)	Reference Ranges (m)	Reference Ranges (a)
LDL-Chol/HDL-Chol	0.716-3.333	0.824-4.768	0.756-4.469
Apo-B100/Apo-A1	0.220-0.926	0.250-1.181	0.232-1.118
Free LDL-Chol/LDL-Chol	0.207-0.208	0.207-0.208	0.207-0.208
Free-HDL-Chol/HDL-Chol	0.151-0.166	0.153-0.176	0.152-0.174

Particle size

Parameter	Reference Ranges (f) nm	Reference Ranges (m) nm	Reference Ranges (a) nm
Chylomicrons	80.370-86.859	80.370-87.235	80.370-87.139
Total VLDL-particle size	34.949-36.199	35.216-37.537	35.053-37.284
VLDL-1 particle size	57.960-63.186	58.262-63.395	57.960-63.309
VLDL-2 particle size	49.446-50.534	50.090-50.589	49.876-50.568
VLDL-3 particle size	42.641-42.862	42.708-42.911	42.666-42.892
VLDL-4 particle size	36.477-36.764	36.539-36.833	36.502-36.801
VLDL-5 particle size	31.579-31.685	31.565-31.746	31.574-31.795
Total LDL-particle size	22.148-22.609	22.104-22.537	22.116-22.571
LDL-1 particle size	24.371-24.455	24.366-24.446	24.369-24.451
LDL-2 particle size	22.938-22.948	22.936-22.945	22.937-22.947
LDL-3 particle size	22.350-22.350	22.350-22.350	22.350-22.350
LDL-4 particle size	21.825-21.834	21.826-21.833	21.826-21.834
LDL-5 particle size	21.130-21.136	21.135-21.137	21.131-21.137
LDL-6 particle size	19.550-19.635	19.530-19.611	19.537-19.625
Total HDL-particle size	8.644-9.186	8.441-9.078	8.487-9.157
HDL-1 particle size	11.422-11.519	11.427-11.613	11.423-11.581
HDL-2 particle size	10.221-10.228	10.214-10.226	10.216-10.227
HDL-3 particle size	9.612-9.637	9.595-9.636	9.601-9.637
HDL-4 particle size	8.172-8.361	8.130-8.315	8.138-8.351

*The reference values are derived from the 95% concentration range of 3242 measured samples from healthy donors (1591 female and 1651 male donors, > 18 years old) whereby the sample handling took place under very controlled conditions. Reference ranges specified by gender, female (f), male (m). If no gender was specified, the gender-neutral reference range for all persons (a) is used.

Particle number

Parameter	Reference Ranges (f) nmol/L	Reference Ranges (m) nmol/L	Reference Ranges (a) nmol/L
Non-HDL particle number	480-1531	469-1664	475-1605
LDL- particle number	336-1198	334-1169	418-1396
IDL- particle number	30-63	28-69	29-66
VLDL- particle number	17-126	23-189	19-170
HDL- particle number	11648-21160	10027-19011	10428-20563
Chylomicrons-particle number			<1
HDL-1 particle number	840-3165	313-2548	425-2969
HDL-2 particle number	581-1922	318-1544	374-1828
HDL-3 particle number	1056-2677	794-2182	858-2541
HDL-4 particle number	8514-13869	8038-13063	8221-13613
LDL-1 particle number	180-477	171-492	175-484
LDL-2 particle number	61-188	60-199	61-192
LDL-3 particle number	21-73	21-78	21-75
LDL-4 particle number	29-114	30-123	29-118
LDL-5 particle number	24-99	25-107	24-104
LDL-6 particle number	104-406	103-447	103-430
VLDL-1 particle number	<3	<6	<4
VLDL-2 particle number	<6	<11	<9
VLDL-3 particle number	2-17	3-30	2-25
VLDL-4 particle number	8-47	9-73	8-64
VLDL-5 particle number	7-54	9-74	8-69

The reference values are derived from the 95% concentration range of 3242 measured samples from healthy donors (1591 female and 1651 male donors, > 18 years old) whereby the sample handling took place under very controlled conditions. Reference ranges specified by gender, female (f), male (m). If no gender was specified, the gender-neutral reference range for all persons (a) is used.

Lipoprotein Main Classes

Parameter	Reference Ranges (f) mg/dL	Reference Ranges (m) mg/dL	Reference Ranges (a) mg/dL
Cholesterol in			
HDL	52.3-123.9	38.6-103.7	41.4-117.8
LDL	70.1-208.5	66.9-221.2	68.5-215.3
IDL	10.5-22.7	10.1-25.0	10.3-23.8
VLDL	4.6-38.4	6.3-63.4	5.1-53.9
Chylomicrons	<1.3	<3.8	<2.7
Esterified Cholesterol in			
HDL	73.5-176.7	53.5-147.0	57.7-167.6
LDL	93.4-277.6	89.2-294.4	91.3-286.7
IDL	14.4-31.1	13.7-34.1	14.0-32.5
VLDL	5.8-48.5	7.9-81.0	6.4-68.4
Chylomicrons	<2	<5.6	<4.0
Phospholipids in			
HDL	73.0-158.6	56.7-103.2	60.9-152.4
LDL	37.7-118.4	36.4-126.4	37.0-122.6
IDL	7.3-15.6	6.9-17.2	7.1-16.4
VLDL	4.2-32.8	5.8-51.8	4.7-45.3
Chylomicrons	<0,1	<0,2	<0.1
Free Cholesterol in			
HDL	8.7-18.8	6.7-16.1	7.2-18.1
LDL	14.5-43.4	13.9-46.0	14.2-44.8
IDL	2.0-4.3	1.9-4.3	1.9-4.5
VLDL	1.1-9.6	1.6-15.3	1.3-13.2
Chylomicrons	<0.2	<0.5	<0.3
Triglycerides in			
HDL	8.1-19.4	6.0-16.2	6.4-18.4
LDL	11.7-36.7	11.2-39.2	11.4-38.02
IDL	13.8-29.8	13.2-32.7	13.5-31.2
VLDL	21.6-176.2	29.3-290.0	23.7-247.4
Chylomicrons	<6.1	<17.7	<1.3
Apo-A1 in			
HDL	147.2-269.3	125.7-240.8	130.4-262.1
Apo-A2 in			
HDL	31.4-54.5	27.9-49.7	29.1-53.1
Apo-B100 in			
LDL	41.6-139.5	40.6-150.3	40.7-145.6
IDL	4.0-8.6	3.8-9.4	3.9-9.0
VLDL	2.0-14.8	2.7-22.0	2.2-19.9
Chylomicrons			<0.1

The reference values are derived from the 95% concentration range of 3242 measured samples from healthy donors (1591 female and 1651 male donors, >18 years old) whereby the sample handling took place under very controlled conditions. Reference ranges specified by gender, female (f), male (m). If no gender was specified, the gender-neutral reference range for all persons (a) is used.

LDL-Subclasses

Parameter	Reference Ranges (f) mg/dL	Reference Ranges (m) mg/dL	Reference Ranges (a) mg/dL
Cholesterol in			
LDL-1	38.3-100.6	36.2-103.9	37.3-102.4
LDL-2	10.4-31.7	10.1-31.7	10.3-32.7
LDL-3	3.2-11.1	3.2-12.0	3.2-11.5
LDL-4	4.1-15.9	4.3-17.2	4.2-16.5
LDL-5	3.0-13.0	2.9-12.7	3.0-12.8
LDL-6	9.8-37.9	9.7-41.7	9.7-40.1
Esterified Cholesterol in			
LDL-1	51.2-134.6	48.5-139.1	49.9-137.0
LDL-2	13.8-42.2	13.5-44.7	13.7-43.1
LDL-3	4.3-14.8	4.4-15.9	4.3-15.2
LDL-4	5.4-21.1	5.7-22.8	5.5-21.8
LDL-5	3.9-16.3	4.0-17.6	4.0-17.1
LDL-6	12.9-50.3	12.8-55.2	12.8-53.1
Phospholipids in			
LDL-1	16.7-43.4	15.8-44.8	16.2-44.2
LDL-2	4.5-13.7	4.4-14.5	4.4-14.0
LDL-3	1.7-6.0	1.7-6.4	1.7-6.1
LDL-4	2.7-10.8	2.9-11.7	2.8-11.1
LDL-5	2.4-10.1	2.5-10.9	2.5-10.6
LDL-6	9.0-35.2	8.9-38.7	8.9-37.1
Free Cholesterol in			
LDL-1	7.8-20.6	7.4-21.3	7.6-20.9
LDL-2	2.2-6.7	2.1-7.1	2.2-6.8
LDL-3	0.7-2.4	0.7-2.5	0.7-2.4
LDL-4	0.9-3.4	0.9-3.7	0.9-3.5
LDL-5	0.6-2.6	0.7-2.8	0.6-2.8
LDL-6	2.1-8.0	2.0-8.8	2.1-8.5
Triglycerides in			
LDL-1	5.1-13.2	4.8-13.6	4.9-13.4
LDL-2	1.3-4.0	1.3-4.3	1.3-4.1
LDL-3	0.5-1.8	0.5-1.9	0.5-1.8
LDL-4	0.8-3.3	0.9-3.6	0.8-3.4
LDL-5	0.8-3.2	0.8-3.4	0.8-3.3
LDL-6	2.9-11.5	2.9-12.6	2.9-12.1

The reference values are derived from the 95% concentration range of 3242 measured samples from healthy donors (1591 female and 1651 male donors, >18 years old) whereby the sample handling took place under very controlled conditions. Reference ranges specified by gender, female (f), male (m). If no gender was specified, the gender-neutral reference range for all persons (a) is used.

HDL-Subclasses

Parameter	Reference Ranges (f) mg/dL	Reference Ranges (m) mg/dL	Reference Ranges (a) mg/dL
Cholesterol in			
HDL-1	12.9-47.7	5.3-38.6	6.8-44.9
HDL-2	5.1-16.9	2.8-13.6	3.3-16.1
HDL-3	6.9-17.6	5.1-14.3	5.6-16.6
HDL-4	24.9-42.4	22.7-39.1	23.3-41.5
Esterified Cholesterol in			
HDL-1	19.2-70.6	7.9-57.2	10.1-66.4
HDL-2	7.3-24.3	4.0-19.5	4.7-23.1
HDL-3	9.7-24.9	7.2-20.3	7.9-23.5
HDL-4	33.7-57.7	30.8-53.0	31.5-56.4
Phospholipids in			
HDL-1	12.5-46.3	5.0-37.5	6.4-43.5
HDL-2	6.2-20.6	3.4-16.5	4.0-19.5
HDL-3	9.6-24.4	7.2-19.9	7.8-23.2
HDL-4	40.8-69.2	37.6-64.0	38.5-67.6
Free Cholesterol in			
HDL-1	1.5-5.7	0.6-4.6	0.8-5.4
HDL-2	0.7-2.4	0.4-2.0	0.5-2.3
HDL-3	1.1-2.8	0.8-2.7	0.9-2.6
HDL-4	4.8-8.1	4.4-7.5	4.6-7.9
Triglycerides in			
HDL-1	2.0-7.4	0.8-6.0	1.1-7.0
HDL-2	0.8-2.6	0.4-2.1	0.5-2.8
HDL-3	1.1-2.8	0.8-2.3	0.9-2.7
HDL-4	3.8-6.6	3.5-6.1	3.6-6.5

*The reference values are derived from the 95% concentration range of 3242 measured samples from healthy donors (1591 female and 1651 male donors, >18 years old) whereby the sample handling took place under very controlled conditions. Reference ranges specified by gender, female (f), male (m). If no gender was specified, the gender-neutral reference range for all persons (a) is used.

VLDL-Subclasses

Parameter	Reference Ranges (f) mg/dL	Reference Ranges (m) mg/dL	Reference Ranges (a) mg/dL
Cholesterol in			
VLDL-1	<3.9	<8.3	<6.4
VLDL-2	0.6-9.2	0.6-9.2	0.4-7.6
VLDL-3	0.8-8.2	1.2-14.0	0.9-11.6
VLDL-4	2.0-12.3	2.5-19.4	2.1-16.8
VLDL-5	1.3-9.7	11.7-13.4	1.4-12.5
Esterified Cholesterol in			
VLDL-1	<5.5	<11.6	<9.0
VLDL-2	0.4-6.7	0.8-12.3	0.6-10.2
VLDL-3	1.0-10.5	1.5-18.0	1.2-14.9
VLDL-4	2.5-15.3	3.1-24.1	2.6-20.8
VLDL-5	1.6-11.5	2.0-15.9	1.7-14.9
Phospholipids in			
VLDL-1	<2.1	<4.4	<3.4
VLDL-2	0.2-3.1	0.4-5.7	0.3-4.7
VLDL-3	0.6-5.7	0.8-9.7	0.6-8.1
VLDL-4	1.9-11.8	2.4-18.5	2.1-16.1
VLDL-5	01.4-10.6	1.9-14.7	1.6-13.7
Free Cholesterol in			
VLDL-1	<0.7	<1.4	<1.1
VLDL-2	0.1-1.0	0.1-1.9	0.1-1.6
VLDL-3	0.2-1.9	0.3-3.3	0.2-2.7
VLDL-4	0.5-3.3	0.7-5.1	0.6-4.4
VLDL-5	0.4-2.8	0.5-3.9	0.4-3.7
Triglycerides in			
VLDL-1	<18.6	<39	<30.3
VLDL-2	1.3-20.8	2.6-38.0	1.8-31.4
VLDL-3	3.6-36.5	5.3-62.4	4.1-51.6
VLDL-4	9.8-60.7	12.2-95.7	10.5-82.6
VLDL-5	5.9-43.4	7.7-60.1	6.5-55.9

*The reference values are derived from the 95% concentration range of 3242 measured samples from healthy donors (1591 female and 1651 male donors, >18 years old) whereby the sample handling took place under very controlled conditions. Reference ranges specified by gender, female (f), male (m). If no gender was specified, the gender-neutral reference range for all persons (a) is used.

References

1. Mora S, Glynn RJ, Ridker PM. High-density lipoprotein cholesterol, size, particle number, and residual vascular risk after potent statin therapy. *Circulation*. 2013;128(11):1189-1197. Doi:10.1161/CIRCULATIONAHA.113.002671.
2. Lawer PR, Akinkuolie AO, Ridker PM, et al. Discordance between circulating atherogenic cholesterol mass and lipoprotein particle concentration in relation to future coronary events in women. *Clin. Chem*: 2017;63(4):870-879. Doi:10.1373/clinchem.2016.264515.
3. Akinkuolie AO, Paynter NP, Padmanabhan L, Mora S: High-density lipoprotein particle subclass heterogeneity and incident coronary heart disease. *Circ Cardiovasc Qual Outcomes*. 2014;Jan;7(1):55-63. doi:10.1161/CIRCOUTCOMES.113.000675.
4. Tehrani DM, Zhao Y, Blaha MJ, et al. Discordance of low-density lipoprotein and high-density lipoprotein cholesterol particle versus cholesterol concentration for the prediction of Cardiovascular disease in patients with metabolic syndrome and diabetes mellitus. *Am J Cardiol*. 2016;117(12):1921-1927. Doi:10.1016/j.amjcard.2016.03.040.
5. Mackey RH, Greenland P, Goff DC, Lloyd-Jones D, Sibley CT, Mora S. High-density lipoprotein cholesterol and particle concentrations, carotid atherosclerosis, and coronary events. *J Am Coll Cardiol*. 2012;60(6):508-516. doi:10.1016/j.jacc.2012.03.060.
6. Otvos JD, Shalaurova I, Freedmann DS, Rosenson RS. Effects of pravastatin treatment on lipoprotein subclasses profiles and particle size in the PLAC-I trial. *Atherosclerosis*. 2002;160:41-48.
7. Khera AV, Demler OV, Adelman SJ, et al. Cholesterol efflux capacity, high-density lipoprotein particle number, and incident cardiovascular events: an analysis from the JUPITER trial (Justification for the use of statins in prevention: An intervention trial evaluating rosuvastatin). *Circulation*. 2017;135(25):2494-2504. Doi:10.1161/CIRCULATIONAHA.116.025678.
8. Cromwell WC, Otvos JD, Keyes MJ, Pencina MJ, Sullivan L, Vasan RS, Wilson PW, D'Agostino RB. LDL Particle Number and Risk of Future Cardiovascular Disease in the Framingham Offspring Study - Implications for LDL Management. *J Clin Lipidol*. 2007 Dec;1(6):583-92. doi: 10.1016/j.jacl.2007.10.001. PMID: 19657464; PMCID: PMC2720529.
9. Sachdeva A, Cannon CP, Deedwania PC, Labresh KA, Smith SC Jr, Dai D, Hernandez A, Fonarow GC. Lipid levels in patients hospitalized with coronary artery disease: an analysis of 136,905 hospitalizations in Get With The Guidelines. *Am Heart J*. 2009 Jan;157(1):111-117.e2. doi: 10.1016/j.ahj.2008.08.010. Epub 2008 Oct 22. PMID: 19081406.
10. Neil J. Stone, Jennifer G. Robinson, Alice H. Lichtenstein, C. Noel Bairey Merz, Conrad B. Blum, Robert H. Eckel, Anne C. Goldberg, David Gordon, Daniel Levy, Donald M. Lloyd-Jones, Patrick McBride, J. Sanford Schwartz, Susan T. Shero, Sidney C. Smith, Karol Watson and Peter W. F. Wilson. 2013 ACC/AHA Guideline on the Treatment of Blood Cholesterol to Reduce Atherosclerotic Cardiovascular Risk in Adults. *Circulation*. Volume 129, Issue 25_suppl_2, 24 June 2014; Pages S1-S45. doi.org/10.1161/01.cir.0000437738.63853.7a
11. Toth PP, Grabner M, Punekar RS, Quimbo RA, Cziráky MJ, Jacobson TA. Cardiovascular risk in patients achieving low-density lipoprotein cholesterol and particle targets. *Atherosclerosis*. 2014 Aug;235(2):585-91. doi: 10.1016/j.atherosclerosis.2014.05.914. Epub 2014 May 22. PMID: 24956532.

12. Otvos JD, Mora S, Shalaurova I, Greenland P, Mackey RH, Goff DC Jr. Clinical implications of discordance between low-density lipoprotein cholesterol and particle number. *J Clin Lipidol*. 2011 Mar-Apr;5(2):105-13. doi: 10.1016/j.jacl.2011.02.001. PMID: 21392724; PMCID: PMC3070150.
13. Melander, O, Shiffman, D, Caulfield, M. et al. Low-Density Lipoprotein Particle Number Is Associated With Cardiovascular Events Among Those Not Classified Into Statin Benefit Groups. *JACC*. 2015 Jun, 65 (23) 2571–2573. doi.org/10.1016/j.jacc.2015.02.077
14. Mora S, Otvos JD, Rosenson RS, Pradhan A, Buring JE, Ridker PM. Lipoprotein particle size and concentration by nuclear magnetic resonance and incident type 2 diabetes in women. *Diabetes*. 2010 May;59(5):1153-60. doi: 10.2337/db09-1114. Epub 2010 Feb 25. PMID: 20185808; PMCID: PMC2857895
15. Mackey RH, Kuller LH, Sutton-Tyrrell K, Evans RW, Holubkov R, Matthews KA. Lipoprotein subclasses and coronary artery calcium in postmenopausal women from the healthy women study. *Am J Cardiol*. 2002 Oct 17;90(8A):71i-76i. doi: 10.1016/s0002-9149(02)02636-x. PMID: 12419483.
16. Mora S, Szklo M, Otvos JD, Greenland P, Psaty BM, Goff DC Jr, O'Leary DH, Saad MF, Tsai MY, Sharrett AR. LDL particle subclasses, LDL particle size, and carotid atherosclerosis in the Multi-Ethnic Study of Atherosclerosis (MESA). *Atherosclerosis*. 2007 May;192(1):211-7. doi: 10.1016/j.atherosclerosis.2006.05.007. PMID: 16765964
17. El Harchaoui K, van der Steeg WA, Stroes ES, Kuivenhoven JA, Otvos JD, Wareham NJ, Hutten BA, Kastelein JJ, Khaw KT, Boekholdt SM. Value of low-density lipoprotein particle number and size as predictors of coronary artery disease in apparently healthy men and women: the EPIC-Norfolk Prospective Population Study. *J Am Coll Cardiol*. 2007 Feb 6;49(5):547-53. doi: 10.1016/j.jacc.2006.09.043. Epub 2007 Jan 22. PMID: 17276177.
18. Song TJ, Cho HJ, Chang Y, Youn M, Shin MJ, Jo I, Heo JH, Kim YJ. Low-density-lipoprotein particle size predicts a poor outcome in patients with atherothrombotic stroke. *J Clin Neurol*. 2015 Jan;11(1):80-6. doi: 10.3988/jcn.2015.11.1.80. Epub 2015 Jan 2. PMID: 25628741; PMCID: PMC4302183
19. El Harchaoui K, Arsenault BJ, Franssen R, Després JP, Hovingh GK, Stroes ES, Otvos JD, Wareham NJ, Kastelein JJ, Khaw KT, Boekholdt SM. High-density lipoprotein particle size and concentration and coronary risk. *Ann Intern Med*. 2009 Jan 20;150(2):84-93. doi: 10.7326/0003-4819-150-2-200901200-00006. PMID: 19153411.
20. Parra ES, Panzoldo NB, Zago VH, Scherrer DZ, Alexandre F, Bakkarat J, Nunes VS, Nakandakare ER, Quintão EC, Nadruz W Jr, de Faria EC, Sposito AC. HDL size is more accurate than HDL cholesterol to predict carotid subclinical atherosclerosis in individuals classified as low cardiovascular risk. *PLoS One*. 2014 Dec 3;9(12):e114212. doi: 10.1371/journal.pone.0114212. PMID: 25470778; PMCID: PMC4254940.
21. Freedman DS, Otvos JD, Jeyarajah EJ, Barboriak JJ, Anderson AJ, Walker JA. Relation of lipoprotein subclasses as measured by proton nuclear magnetic resonance spectroscopy to coronary artery disease. *Arterioscler Thromb Vasc Biol*. 1998 Jul;18(7):1046-53. doi: 10.1161/01.atv.18.7.1046. PMID: 9672064.